

## Product Data Sheet „motion sprint 600-SL med“



|                        |  |
|------------------------|--|
| name of machine:       | motion sprint 600-SL med   |
| manufacturer:          | Sprintex Trainingsgeräte GmbH  |
| part number:           | F-MED-ML-610   |
| monitor:               | E900- monitor with casing, display and overlay   |
| power range:           | 0-17 km/h in 0,1 km/h steps, incline 0-15 %  |
| motor system:          | brushless 1100W DC-Motor, tooth belt controlled lamella system   |
| applied standards:     | DIN EN 20957-1, DIN EN 957-6, DIN EN 20957-1, DIN EN 957-6, DIN EN 60601-1, DIN EN 60601-1-2/Safety Class Power supply: 230V/50-60 Hz/EEC, separate fuse required  |
| classification:        | 93/42 EEC class IIa  |
| usage class:           | S, I following DIN EN ISO 20957-1  |
| accuracy class:        | A (high accuracy) following DIN EN ISO 20957-1   |
| ambient conditions:    | humidity <80 %; ambient temperature 0-45°C   |
| display (resolutions): | LCD and Dot Matrix Display with following parameters: program number (1-10); time (00:00); pulse (optional); watt (1-999 in 1 watt-steps); km/h (0,1-99,9); km (in meter 1-999), after that in 10m-steps 1,00-99,99) ; k-cal (0,0-999,9); level (1-21), incline (0-15) |
| pulse measurement:     | coded Polar® receiver, earpulse optional   |
| interfaces:            | standard: RS232-interface, chipcard system<br>optional: galvanically isolated RS232-interface, RFID, NFC, wireless network. Please consult on individual settings/needs in time.   |
| programs:              | standard: Quick, Watt, Pulse<br>optional: Watt sensitive, Therapy, Interval, Calorie, Hill, Random, Distance   |
| maximum user weight:   | 150 kg   |
| measurements:          | 173 x 76 x 145 cm (l/w/h)  |
| net weight:            | 170 kg   |
| gross weight:          | 203-275 kg (depending packaging and mode of transport)   |

Weight and dimensions may differ depending on accessories.

Frame colour depending on individual specifications.

Technical changes reserved.

Warning: Only authorized service staff is allowed to open the unit.

If you need further information please refer to user manual.

(Issued: May 2021)